



# SFC3A

The SFC3A is easily recognizable by its classic design. Its large gripping area makes it an easy to handle microscope throughout the classroom and even outdoors. This series is designed with the elementary school and hobby users in mind.

<b>Model</b>	SFC-3A
<b>Optical System</b>	Finite Optical System
<b>Observation tube</b>	Monocular
<b>Inclination</b>	Adjustable inclination
<b>Eyepiece(s)</b>	WF10X/14mm with pointer, lockable
<b>Nosepiece</b>	Triple revolving nosepiece with click stops
<b>Objective classification</b>	Achromatic, DIN
<b>Objective information</b>	4X/0.10, 10X/0.25, 40X/0.65/S
<b>Stage</b>	Plain stage with sample clips
<b>Stage size</b>	88X82mm
<b>Condenser</b>	Rotary disc diaphragm
<b>Focus mechanism</b>	Coarse focusing system with tension adjustment
<b>Illumination</b>	Plano-concave mirror and fork
<b>Accessories included</b>	Dust cover
<b>Length X Width X Height</b>	158X111X290
<b>Weight KG</b>	1,3KG